















Structure Search

Document Search

Custom Synthesis / Bulk Chemicals V References V How to Order V About TCI V

FAQ ~

New tool "TCI SpectraViewer" | [TCIPracticalExample] Methyl Esterification of Carboxylic Acid Using...

## CAS RN: Product Number: T3983

Tissue-Clearing Reagent CUBIC-R+(N) [for Animals]



Size	Unit Price	Saitama (Kawaguchi)	Hyogo (Amagasaki)	Stock In Other WH	More Stock Info ②
25ML	¥4,000	17	Contact Us	Contact Us	Check Estimated Ship Date
100ML	¥11,000	3	34	Contact Us	Check Estimated Ship Date
500ML	¥40,000	13	Contact Us	Contact Us	Check Estimated Ship Date

Request Bulk Quote

#### Note:

A tissue clearing method "CUBIC" has been developed by Prof. Hiroki R. Ueda and coworkers at The University of Tokyo / RIKEN. The CUBIC technique enables cyclopedic imaging at a single-cell resolution following whole-body and whole-organ clearing. CUBIC reagents

<sup>\*</sup> Check with "Check Estimated Ship Date" for details. The shipping date may vary depending on the quantity. If it is displayed as 4 weeks or more, it may be shortened depending on the order.

<sup>\*</sup> For the case that without specific estimated shipping data, e.g.: If displayed as "will be shipped after 8 working days", the products will be shipped 8 working days later after ordering.

<sup>\*</sup> Available Stock: Prompt shipment (for products) in Saitama/Hyogo warehouse. Stock in other WH: 2-3 Business Days

<sup>\*</sup> Please contact us if you need further information.

<sup>\*</sup> To send your quote request for bulk quantities, please click on the "Request Bulk Quote" button. Please note that we cannot offer bulk quantities for some products.

<sup>\*</sup>TCI frequently reviews storage conditions to optimize them. Please note that the latest information on the storage temperature for the products is described on our website.

which can be used in the tissue clearing method are provided by TCI.

\*CUBIC reagents are commercialized under invention licenses by RIKEN, Japan.

Two products of CUBIC-L (Product Number: T3740), and CUBIC-R+(N) (Product Number: T3983) or CUBIC-R+(M) (Product Number: T3741) are required for tissue-clearing.

The difference between CUBIC-R+(N) (Product No. T3983) and CUBIC-R+(M) (Product No. T3741):

CUBIC-R+(N) (Product No. T3983) is inexpensive and easier to handle because it raises less precipitation. The fluorescence signal may decay, but the fluorescence signals of samples in CUBIC-R+(N) can be observed for several days after immersion.

CUBIC-R+(M) (Product No. T3741) is superior in retaining the fluorescence signal. However, at low temperatures such as in winter, it may precipitate. In that case, it can be resolved by placing the sample at 37°C for a few days.

For these reasons, it is recommended to try CUBIC-R+(N) (Product No. T3983) first and then use CUBIC-R+(M) (Product No. T3741) if fluorescence signal cannot be found.

Please see product category page or a technical guidebook (PDF) of how to use. And an insertion (PDF) is also available.

General Information Specifications & Properties Safety & Regulations Applications

## Articles/Brochures




CUBIC - Animal Tissue-Clearing Reagents - Technical Guidebook



Animal Tissue-Clearing Reagent: CUBIC

### **TCIMAIL**

Documents (Note: Some products will not have analytical charts available.)



# **Product Categories**

#### Life Science

Biochemicals and Reagents > Tissue-Clearing Reagents > Animal Tissue-Clearing Reagents

Cell Biology > Cell and Tissue Imaging > Animal Tissue-Clearing Reagents

**How To Order** 

Company Profile

**About TCI** 

How to Order Our Distributor TCI Philosophy

CSR Request Bulk Quote

Custom Synthesis Request Our Story

**Product Document Search** 

Ordering FAQ Locations Sign up for TCI eNews & Promotions

**Subscribe** 

Subscribe above to receive our eNews, featuring scientific research, new products, and exclusive promotions. You can easily







unsubscribe at any time. Sign up today!





Contact Us

©Tokyo Chemical Industry Co., Ltd.

**Privacy Policy** 

**Terms & Conditions**